



Health-Shine Technology Co.,Ltd

Health -Shine Technology Co.,Ltd (HST) was founded in the year of 2008,which is a professional manufacturer and supplier of Spo2 sensor and adapter , patient monitor ECG cable and leadwires , Temperature probe , NIBP cuff and hose , Electrodes and various medical instruments. HST products are complied with ISO 13485-2003 and latest medical international standards.

With guidance of Quality First and Customer Supremacy,we have established stable relationship with clients from all of the world.Our products have good demand in South of Asia,Middle East ,Europe and America. We supply the best quality and competitive price to our partner to success in the marketing.

HST has following advantages to help you to be successful

?.Professional R&D technical team support

?.Qualityied ,Stable Product and Competitive Price

?.Timely Delivery

Welcome to visit Health-Shine Technology Co.,Ltd.

Contact:

Mr.Andy



Tel:+86-755-27633663

Fax:+86-755-27633668

Email:andy@hstmed.com

MSN:hstmed@hotmail.com

skype:hstmed

Original Accessories



[Masimo spo2 sensor LNOP DCI](#)

[Nellcor DS-100A adult finger spo2](#)

[Tyco Nellcor spo2 sensor original sensor](#)

IBP Series



[Philips-Utah IBP cable](#)



[Datex- Abbott IBP Cable](#)



[GE-Appott IBP cable](#)



[Spacelabs-Appott IBP cable](#)



[Edwaeds IBP Pressure Trandusers](#)



[Abbott IBP Pressure Trandusers](#)

NIBP Series



[NIBP connector](#)



[Disposable NIBP Cuff](#)



[Reusable Cuff with bladder](#)



[Reusable Cuff without bladder](#)



[Datex NIBP hose](#)



[HP/Philips NIBP Air Hose](#)



[GE NIBP Hose](#)

Temp.Series



[YSI700 Skin Temp. Probe](#)



[YSI400 Esophageal-Rectal
Temp.Probe](#)

EKG Cable



[\(SH\)Nihon Kohden EKG cable with 10 leads](#)



[KENZ EKG cable with 10 leadwire](#)



[GE-MAC EKG cable with 10 leadwire](#)



[HP/Philips EKG Cable](#)



[Fukuda EKG cable with 10 leadwire](#)



[Schiller EKG cable with 10 leads](#)

ECG Cable



[Multi-Functional ECG Limb Clamps](#)



[Adult ECG Sunction Electrode](#)



[HP/PHILIPS ECG 3-Lead Trunk Cable M1510A](#)



[GE ECG cable with 3 leadwires](#)



[GE ECG cable with 3 leadwires](#)



[HP/Philips ECG Cable 12Pin](#)



[HP/Philips ECG Cable 8Pin](#)



[AAMI ECG cable 6Pin](#)



[HP/Philips ECG leadwires AA](#)



[Siemens ECG trunk cable](#)



[Datex ECG cable](#)

SpO2 Sensor



[Nellcor OxiMaxSensor](#)



[Nellcor spo2 adapter cable](#)



[HP/Philip spo2 adapter cable](#)



[Nellcor extension cable DEC-4](#)



[HP/Philip spo2 adapter cable](#)



[Masimo spo2 adapter cable](#)



[Nellcor Pediatric spo2 sensor](#)



[Masimo spo2 sensor](#)



[GE Marquette Adult finger spo2 sensor](#)



[Masimo adult spo2 sensor](#)



[Nihon Kodhen Adult softtip spo2 sensor](#)



[Nellcor Disposable Neonate spo2 sensor](#)



[Nellcor Disposable Adult spo2 sensor](#)



[Nellcor Disposable infant spo2 sensor](#)



[Masimo disposable Neonate spo2 sensor](#)

Contact Information

Company Name: Health-Shine Technology Co.,Ltd

Contact Person: Andy Shi

Street Address: 5F,139#, Ruyi Road,Shiyan Town,Bao'an District,

City: Shenzhen

Country/Region:

Zip: 518000

Telephone: 86-755-27633663

Mobile Phone: +8613823751550

Fax: 86-755-27633668

Email: andy@hstmed.com

Company Site: <http://hstmed.medincn.com>